

## Adequate jacket size for Sprague Dawley rats :

		Body weight range of <b>male</b> rat in grams*										Dressing tube	
		180-230	210-290	250-330	280-370	310-400	330-420	350-440	370-470	390-490	400-540		420-650
Jacket size	0	> 220	250-270										Small
	1		250-270		350-370								Small
	2				350-370				450-470				Small / Medium
	3								450-470				Medium / Large
		7	8	9	10	11	12	13	14	15	16	17 +	
<b>Age in weeks**</b>													

		Body weight range of <b>female</b> rat in grams*									Dressing tube	
		190-230	200-250	210-270	220-280	230-290	235-295	240-300	250-310	260-320		270-400
Jacket size	0	> 220						290-300				Small
	1							300-310				Small / Medium
		9	10	11	12	13	14	15	16	17	18+	
<b>Age in weeks**</b>												

\* based on the growth chart of the Outbred CD® (Sprague Dawley) IGS rat from Charles river (breeding place: Germany, UK and Italy).

\*\* The age (week) is an additional approximate indication. The jackets are to be chosen according to the weight of the animals.

## Adequate jacket size for Wistar rats :

		Body weight range of <b>male</b> rat in grams*										Dressing tube	
		190-240	220-270	250-310	270-330	290-370	300-390	320-410	330-430	450-470			
Jacket size	0	> 220		250-270									Small
	1			250-270			350-370						Small
	2						350-370			450-470			Small / Medium
	3									450-470			Medium / Large
		7	8	9	10	11	12	13	14	15	16	17 +	
<b>Age in weeks**</b>													

		Body weight range of <b>female</b> rat in grams*									Dressing tube	
		170-230	190-240	200-260	210-260	210-270						
Jacket size	0	> 220					290-300					Small
	1						300-310					Small / Medium
		10	11	12	13	14	15	16	17	18+		
<b>Age in weeks**</b>												

\* based on the growth chart of the Outbred CD® (Wistar) IGS rat from Charles river (breeding place: Germany, UK and Italy).

\*\* The age (week) is an additional approximate indication. The jackets are to be chosen according to the weight of the animals.