

SurgiSuite

Multi-functional
surgical platform

ANIMA LAB
animal facility and laboratory equipment • animal research models



Chemical Resistant Base

*Stainless steel surgical surface
and polycarbonate overlay for
added protection*

Far Infrared Animal Warming Pad (optional)

*Recessed in center of
platform for a level
surgical field*

Bright LED Lighting (optional)

*Connects to multi-port
USB power strip*

SurgiSuite shown with optional add-ons

kentscientific.com

Order your SurgiSuite:
www.animalab.eu | info@animalab.eu



Kent Scientific
CORPORATION

SurgiSuite includes:

- Easy to clean, chemical-resistant surface
- Magnetic anesthesia mask/intubation tubing stabilizer, (1-small and 1-large)
- Magnetic limb positioners (4)
- Elastomer wire positioner (1-meter)
- Detachable procedure bowl
- Disposable platform covers (5-pack)

1-Year 100% Parts & Labor Warranty

Visit our website for a complete selection of SurgiSuite accessories.

Specifications:

	SurgiSuite - Standard	SurgiSuite - Large
Base dimensions:	12.90 in x 16.80 in (32.77 cm x 42.67 cm)	15.75 in x 17.80 in (40.01 cm x 45.21 cm)
Stainless steel platform:	11.90 in x 7.90 in (30.23 cm x 20.07 cm)	12.75 in x 10.75 in (32.39 cm x 27.31 cm)
Warming surface:	6.6 in x 5 in (16.51cm x 12.7cm)	8.75 in x 6.75 in (22.23 cm x 17.15 cm)

Optional Add-Ons:



Monitor and regulate animal temperature

The PhysioSuite® with RightTemp® temperature monitor and homeothermic control module easily connects to the SurgiSuite.

Ask about integrating pulse oximetry.



Customize your lighting solution

LED lighting plugs directly into the included USB power strip.